

Detect Lysosomal Iron

DOJINDO

with Precision

► Lyso-FerroRed

This product was developed under the guidance of Professor Tasuku Hirayama, Gifu Pharmaceutical University.




How Crucial Is Lysosomal Iron in Ferroptosis?

Ferroptosis, an iron-dependent form of cell death driven by lipid peroxidation, is implicated in cancer, neurodegenerative, and age-related diseases. While early studies focused on total cellular iron, increasing attention has shifted to its subcellular distribution, particularly within lysosomes. As major sites of iron storage, lysosomes may provide a critical trigger for ferroptosis*, highlighting the importance of technologies that enable selective detection and evaluation of lysosomal iron.


*Saimoto, Y, et al., *Nat Commun* **2025** 16, 3554.

Iron Redistribution and Ferroptosis Upon Stimulation




Product Feature


Lyso-FerroRed selectively reacts with Fe^{2+} in lysosomes, producing strong fluorescence and clearly revealing lysosomal localization. Co-staining with a lysosomal dye (LysoPrime Deep Red, Cat#L264) showed overlapping signals, confirming its lysosome-specific staining.



*The merged image is shown in pseudo-color for clarity.



Workflow



Feature Comparison

Flexible selection according to your assay and instrument

	Lyso-FerroRed	FerroOrange	Mito-FerroGreen
Localization	Lysosome	Intracellular*	Mitochondria
Fluorescence	$\lambda_{\text{ex}}=551 \text{ nm}$, $\lambda_{\text{em}}=571 \text{ nm}$	$\lambda_{\text{ex}}=543 \text{ nm}$, $\lambda_{\text{em}}=580 \text{ nm}$	$\lambda_{\text{ex}}=505 \text{ nm}$, $\lambda_{\text{em}}=535 \text{ nm}$
Instruments	Fluorescence microscope, Plate reader, FCM	Fluorescence microscope, Plate reader	Fluorescence microscope
Sample	Live Cells	Live Cells	Live Cells


*For information on the staining localization, please visit our website.

NEXT » Featured Experiment

Featured Experiment


Fluorescence Microscopy Imaging

Experiments were performed in HeLa cells, and changes in lysosomal Fe²⁺ were monitored with Lyso-FerroRed in the presence or absence of added iron (ferrous ammonium sulfate, final concentration: 200 μM) and the iron chelator Bpy (2,2'-bipyridine, final concentration: 1 mM).



Effects of Ferroptosis Inducer and Lysosomal Inhibitor Co-treatment

Ferroptosis-resistant A549 cells were treated for 24 h with the ferroptosis inducer RSL3, alone or in combination with the lysosomal inhibitor chloroquine (CQ). RSL3 alone caused minimal effects on viability or lysosomal Fe²⁺, though partial lysosomal accumulation was observed. In contrast, co-treatment with CQ increased lysosomal Fe²⁺ and lysosomal size while reducing viability, suggesting a role of lysosomal Fe²⁺ in ferroptosis promotion.



Product name	Unit	Code
Lyso-FerroRed	35 nmol	L270-10

Typical Usage: 17 x 35 mm dishes per 35 nmol

Product name	Unit	Code
FerroOrange	1 tube	F374-10
	3 tubes	F374-12
Mito-FerroGreen	50 μg x 2	M489-10
LysoPrime Deep Red - High Specificity and pH Resistance	1 tube	L264-10
	3 tubes	L264-12

Typical Usage: FerroOrange: 17 x 35 mm dishes per 1 tube; Mito-FerroGreen: 5 x 35 mm dishes per 50 μg; LysoPrime Deep Red: 10 x 35 mm dishes per 1 tube

Distributed by

Contact to

Dojindo Europe GmbH

Leopoldstr. 254

80807 Munich, Germany

Email: info@dojindo.eu.com

Manufactured by

DOJINDO

in Japan