



PRODUCT INFORMATION

96 Well Filter Plate (960ul each Well)

Product name: 96 Well Filter Plate
Catalogue Number: SD5006

Description:

The product is a colorless 96 Well Filter Plate with white filter membrane and blue O-ring. Our 96 Well Filter Plate offers a convenient way to remove proteins from genomic and plasmid DNA in a large scale format. They are able to remove viscous proteins prior to DNA binding step.

Technical Specifications:

	96 Well Filter Plate
Color	Natural/clear
Plastic materials	Polypropylene
Filter membrane	Yes
Material of Membrane	Polypropylene
Pore size of the filter	10 μ m
Thickness of the filter	Around 2.5 mm
Color of filter membrane	White
O-ring	Blue
Dimensions (W x D x H)	128.2mm x 84.9mm x 37.5mm
Weight	105g to 120g
Package	12/bag

Storage

Keep product tightly closed in a dry, cool and well-ventilated place.

Intended for basic research purpose only.