

PRODUCT INFORMATION

0.1-10µl Filter Tips

Description:

0.1-10UL filter tip.

Pipette tips with filter inserts avoid contamination with liquids and aerosols during the pipetting. Filter tips minimize the danger of unwanted cross-contamination.

Features

- Filter tip. Universal fit: Eppendorf, Gilson and most popular pipettes.
- 0.1-10UL
- Safely and accurately pipette
- Autoclavable
- Storage RT



Item No.	Description	Package
FT215	Natural/filter/bag	100/bagx10
FT215-RNS	Natural/filter/rack/sterile Certificated RNase, DNase-Free.	96/rack
FT215-RNSB	Natural/filter/rack/sterile Certificated RNase, DNase-Free.	96/rackx10

Technical Specifications

Item No.	FT215	FT215-RNS	FT215-RNSB
Package	100/Bagx10	96/Rack	96/Rackx10
Volume	0.1-10µl	0.1-10µl	0.1-10µl
Graduations	NO	NO	NO
Color	Natural	Natural	Natural
Plastic materials	Polypropylene	Polypropylene	Polypropylene
Filter	YES	YES	YES
Sterile	NO	YES	YES
DNase,RNase Free	NO	YES	YES
Withstand Temperatures	-80°C~121°C	-80°C~121°C	-80°C~121°C

Comparison of Related Tips



QF 24 Rev 1